


# Keep It Clean

## Overflow Safety Traps & Suction Filters Help Prevent Contamination

Overflow Safety Traps and Suction Filters are simple, inexpensive devices that protect your vacuum regulators and vacuum systems from accidental ingress of contaminants.

Filters and Traps help reduce maintenance costs by decreasing the need to clean and disinfect your vacuum regulators because of contamination.

As a preventative measure, keeping them clean is less expensive than cleaning, or replacing them after the damage has been done.



*Interior View of Push-to-Set™ Intermittent Suction Unit*


**Don't let your Vacuum Regulators get like this!**



+



=



**Use an Overflow Safety Trap & Suction Filter!**

# Your Final Lines of Defense

## Overflow Safety Traps



### Overflow Safety Traps

The overflow safety trap guards against accidental overflow and protects regulators, wall outlets and central pipeline systems from contamination and blockage.

Part Number	Part Description
6700-0365-901	Locking Gland
6702-0365-901	DISS Handtight
6703-0365-901	DISS Wing Nut & Gland

## Sure-Traps™ & Filters



**Sure-Trap™**  
There is no cleaning required—simply discard when the trap is contaminated.



**Hydrophobic Filter**  
The disposable high flow hydrophobic filter is designed to reduce contamination of the vacuum regulator by preventing the penetration of fluids and acts as a shut off in the presence of flooding.

Part Number	Part Description
7725-0365-912	12 traps
7725-0365-925	25 traps
7725-0365-950	50 traps

Part Number	Part Description
6730-0570-800	3 filters
6730-0571-800	10 filters
6730-0572-800	50 filters

© 2021 Ohio Medical LLC. This document contains information that is proprietary and confidential to Ohio Medical LLC. Use of this information is under license from Ohio Medical LLC. Any use other than that authorized by Ohio Medical LLC is prohibited. Ohio Medical, the Ohio Medical logo, Amvex, and the Amvex logo are registered trademarks of Ohio Medical LLC.



**Ohio Medical LLC**  
1111 Lakeside Drive, Gurnee, IL 60031 USA  
866.549.6446 | [ohiomedical.com](http://ohiomedical.com)  
255529 (Rev. 2) 07/2021